

## Agfa Cesium DR Panel

\$34,995

### **Includes:**

- Agfa cesium 14x17 or 17x17 Flat Panel Detector \*Your Choice\*
- Desktop Computer
- Musica 3+ PACS
- Three Year Warranty
- Shipping
  
- Picture perfect first-up display
- Large 17×17, 14×17, digital sensor size options
- Large touch control screen
- On-board positioning guide
- Instant email function
- Customizable user interface
- Imaging Technology
- Musica 3+
- AED (Automatic Exposure Detection) Technology
- 10,000,000,000 instructions per image
- Veterinary specific processing
- Advanced Technology
- Thick CSI (Cesium) Low Dose Scintillator
- Waterproof

### **General Details:**


AGFA Healthcare award-winning MUSICA® imaging technology has just released their all-new MUSICA 3® image processing. Graduating from the older CR processing, the New MUSICA 3 was designed from the ground up specifically for AGFA digital imaging panels for the most consistent diagnostic images available world-wide.

AGFA MUSICA 3+ image processing for companion animals provides exquisite detail throughout the image with specialty algorithms to meet the most demanding clinical needs, processing 10,000,000,000 instructions per image and 1,000 instructions per pixel!


By focusing down to the pixel level, the New MUSICA 3+ image processing, along with the AGFA digital imaging sensors, allows for the most balanced and consistent image quality with every exposure. MUSICA 3+ processing technology provides equal detail in both bone and soft tissue images—all while incorporating the ALARA (As Low As Reasonably Achievable) principles of radiation safety.



**Sample Images**







71