

JPI DirectVet 110V DR System

\$29,995

Includes:

- Heavy Duty Table
- 120kv Generator
- CareRay cesium 14x17 or 17x17 Flat Panel Detector ***Your Choice***
- Laptop or Desktop Computer ***Your Choice***
- ExamVue PACS
- **Five Year Warranty**
- Shipping
- Training

General Details:

- 120 KVP MAXIMUM OUTPUT
- THE HEAVY-DUTY TABLE HOLDS UP TO 300 LBS.
- THE GENERATOR AND DETECTOR MOVE TOGETHER FOR POSITIONING
- SEND DICOM OR JPEG IMAGES
- EMAIL IMAGES FOR REVIEW
- USER FRIENDLY
- SIMPLE AND FAST WORKFLOW
- INTEGRATES WITH YOUR PRACTICE MANAGEMENT SOFTWARE - SUCH AS CORNERSTONE, AVIMARK, IMPROMED, NEO AND MANY MORE!

ExamVue Duo Software Is Specially Designed for the Veterinary Market

- Optimized for Simple and Fast Workflow
- Includes Tools Valuable for Both Large Companies and Small Offices

- Position Instructions and X-ray Techniques for Dogs, Cats, Horses, and Other Animal Exams
- Designed for Touchscreen and Desktop
- Integration with Leading generator manufacturers such as Cubex, Poskom, EcoRay, and Min-X
- ExamVue Duo Software produces clear, high-quality images.
- ExamVue Duo Software comes standard with veterinary diagnostic tools and can be expanded even further with the JPI Cloud PACS.
- ExamVue Duo Software seamlessly integrates into any professional practice

Tethered Flat Panel Detector


14x17 or 17x17 tethered flat panel Cesium panel

Pixel Pitch: 140 μm

Full-field AED for broad system compatibility without the need for sync

Direct-Deposit Cesium

Durable carbon fiber and aluminum alloy housing



Heavy duty table

Holds Over 300 Pounds

110 volt power

The Generator arm slides side to side with the table top

32'D x 54" W

40 inch SID



CubeX 28 X-Ray Generator

2.8 kW

100V-120V / 220V-240V

Single / 50 / 60 Hz

40kV - 120 kV

0.4mAs - 100 mAs



Sample Images







