

Canon 55g series DR Panel

Small Animal DR System

\$14,500

Includes:

- Plate size 14x17 with ideal for the small animal veterinary environment
- Durable HD system computer
- Canon Software included
- High resolution convenient screen interface
- Seamless DICOM® integration to any DICOM conforming PACS system
- X-ray generator sync included
- Enhanced sensitivity for higher quality thick body part imaging
- User interface and styles can be customized to specific needs to help expedite exams.
- Auto-sequence program allows for rapid sequential studies.
- Images appear in 1-2 seconds
- Specific defined animal protocols
- Complete on-image series descriptions and positional markers placed in the proper standard positions, making for accurately identified and labeled studies.

- Use with memory sticks, printers, external drives
- Simple and user-friendly setup









