

Canon 60g series DR Panel

Portable Equine DR System

\$17,650

Includes:

- Plate size 60Gv 9" x 11 with Toughcase packing for reliability and rugged design that is ideal for the veterinary environment
 - Tunnel
 - Software
 - HD Laptop
 - Case
 - Online Setup and Training
 - Shipping Included
 - Optional Cart with folding, adjustable height for secure operation and mobility
 - Durable system case with padded handles and convenient wheels
 - High resolution computer with convenient touch screen or 2 Megapixel non-touch screen interface
 - Seamless DICOM® integration to any DICOM conforming PACS system
 - Optional wireless X-ray generator sync
 - Enhanced sensitivity for higher quality thick body part imaging (stifles, shoulders and necks)
-
- User interface and styles can be customized to specific needs to help expedite exams.
 - Auto-sequence program allows for rapid sequential studies.
 - Images appear in 3 seconds
 - Equine specific defined protocols including well known pre-purchase protocols

- Complete on-image series descriptions and positional markers placed in the proper standard positions, making for accurately identified and labeled studies.
- The system features the sunlight- viewable 1,000 nit touchscreen making it ideal for client presentations, office use and in-the-field use.
- Use with memory sticks, printers, external drives
- Wireless battery powered or hard- wired sync cable 100KV/20mA and 100KV/40mA generators available. Dual laser pointers for accurate SID.



