

Ultrasound Info Sheet  
DP-10

\$3,995

Features:

- 12.1 inch high definition LED
- Full Screen design for easy cleaning
- 30 degree tilting angle adjustable monitor
- User-friendly control panel with backlight
- Light and compact design for extreme portability
- 1.5 hour scanning with rechargeable battery
- Tissue Harmonic Imaging: Significantly enhances contrast resolution and improves image quality
- Tissue Specific Imaging: Optimizes the image quality based on properties of the tissue being scanned. Four options are available including general, muscle, fluid and fat
- iClear™(Speckle Suppression Imaging): Sharper and continuous edges, Smooth and uniform tissues, Cleaner 'clutter-free areas'
- ExFOV: Extended view of the anatomical structure on all convex and linear transducers
- iStorage: Directly transfer images and reports to PC via network cable
- iZoom: Gain instant full screen view on the click of a single key
- iTouch: Obtain instant automatic image optimization in B, Color and PW with one single click
- iStation: Mindray's unique Patient Information Management System allows you to integrate, review, archive and retrieve patient data effectively

