

Ultrasound Info Sheet  
DP-30

\$5,495

Features:

- Mindray DP-30 Veterinary Portable Black and White Ultrasound System
- 12.1" High resolution LED with tilt functionality for better viewing
- One key function for image and cine storage to local disk or USB drive
- Tissue Harmonic Imaging enhancing contrast resolution
- One key quick image optimization with IP (Image Processing)
- 8- segment TGC assuring accurate image adjustment
- Compact size weighing around 12 pounds
- Ready to scan within 17 seconds
- user defined keys—programmable personalized operations
- 1.5 hours continuous scanning
- iStation integrated patient data storage, archiving, review, and retrieval
- iStorage: Directly transfer images and reports to PC via network cable
- iZoom: Gain instant full screen view on the click of a single key
- iTouch: Obtain instant automatic image optimization in B, Color and PW with one single click
- iStation: Mindray's unique Patient Information Management System allows you to integrate, review, archive and retrieve patient data effectively

