

Portable Generators

Large or Small Animal Digital xray System

EcoRay

Plug-in and Hybrid (Battery and A/C)

Your Choice of any EcoRay

\$7,695

BLUE Series: Hybrid Plug-In and Battery Powered Veterinary xray imaging generators, All in One System

ORANGE Series: 110v Veterinary Imaging Generators

1. Control unit
2. High voltage generator
3. X-ray tube
4. Collimator
5. Exposure switch
6. Power cord
7. Carrying case with user manual





Cubex:

40-120 kv

40mA

110v only

26.5lbs

Soft-touch controls with digital displays

Calibrated cassette size indicator dials

Inverted control panel for table use available




Two-stage dynamic, auto line compensation

Easy integration with a digital system

Dual integrated laser pointer



CUBEX

Model		A/C Powered		HYBRID: Battery or A/C Powered	
		9020HF	1040HF	Ultra 90BT	Ultra 1025BT
					
Generator Type		High Frequency Inverter			
HV Generator	Max. Output Power	1.35kW (90kV@15mA)	2.4kW (80kV@30mA)	1.35kW (90kV@15mA)	1.6kW (80kV@20mA)
	Max. kV / mA	90kV / 20mA	100kV / 40mA	90kV / 20mA	100kV / 25mA
X-ray Tube	Focal Spot	1.2 X 1.2mm	1.2 X 1.2mm	0.8 X 0.8mm	1.2 X 1.2mm
Input Power		A/C Power: 110~120V / 60Hz 220~230kV / 50Hz	A/C Power: 110~120V / 60Hz 220~230kV / 50Hz	Battery Power: 60VDC Rechargeable Lithium-Ion A/C Power: 110~120V / 60Hz 220~230kV / 50Hz	Battery Power: 60VDC Rechargeable Lithium-Ion A/C Power: 110~120V / 60Hz 220~230kV / 50Hz
Collimator	Lamp	12V 50W, Halogen Lamp		LED Lamp	
	Dual Laser	Dual laser spot to 66cm target away			
Weight (Approximate)		15 LBS	28 LBS	16 LBS	22 LBS
User Interface	Display	7 Segment LED X2 or 3 digits of kV, mAs Ready, X-ray, Wait, Error LED			
	Operation Switch	<ol style="list-style-type: none"> 1. kV & mAs Up-Down 2. Dual Laser & Lamp ON/OFF 3. Display Reverse 4. APR Memory Storage 5. 8 Preset Memory (90BT and 1025BT have 7 Preset Memory) 			